

Rocket® IPC™ Catheter and Drainage Equipment



Pleural Effusions and Ascites
Rocket® IPC™ Indwelling Pleural & Peritoneal Catheter is indicated for intermittent, long term drainage of symptomatic, recurrent:

- Pleural effusions including malignant pleural effusions.
- Peritoneal ascites, including malignant ascites and other ascites that do not respond to medical management of underlying disease.
- Also the palliation of symptoms related to recurrent malignant ascites.

Rocket® IPC offers an alternative long term method of managing pleural and peritoneal ascites offering significant advantages:

- Early discharge
- Minimised repeat admission.
- Improved quality of life
- Improved control over discomfort and breathlessness
- Reduced cost

DESCRIPTION	ORDER CODE
IPC™ Catheter and Mini Insertion Kit	R54400-MINI
IPC™ DRAINAGE EQUIPMENT Each pack contains a sterile pre-evacuated Drainage Bottle and replacement IPC Cap. All items are single use only.	
IPC™ 1000ml Pre-evacuated Bottle & Cap only	R54411-00-00
IPC™ 500ml Pre-evacuated Bottle & Cap only	R54410
IPC™ ACCESSORIES All items are single use only.	
IPC™ Drainage Connection Line	R54410-00-DL
IPC™ Replacement Cap	R54410-00-CP

- 1 x 16Fg x 40cm Silicone indwelling catheter
- 1 x Drainage line
- 1 x Valve cap
- 1 x Tunneling forceps
- 1 x Foam catheter pad
- 1 x Clear-film dressing
- 1 x 18G Guidewire insertion needle
- 1 x J-Tip Guidewire
- 1 x 16Fr Tearaway sheath/dilator

